




Practi-TRAINER[®]
Universal AED Trainer

WL-120

Quick Start Guide


AED Trainer




Overview

The Practi-Trainer AED Trainer Package

Includes: trainer unit, training pads, wire connector, quick-start guide, and storage case.



Quick Operation Guide



Step 1: Load 3 × AA batteries into the compartment; press On/Off button to power up. Select a scenario on the unit.


Step 2: Follow the voice prompts to place the pads on the manikin and plug in the wire connector.

Step 3: Press the Shock button when advised, then begin CPR as directed.


Descripción general

Paquete del entrenador DEA Practi-Trainer

Incluye: unidad de entrenamiento, parches de entrenamiento, conector de cable, guía de inicio rápido y estuche de almacenamiento.



Guía rápida de funcionamiento



Paso 1: Coloque 3 pilas AA en el compartimento; pulse el botón On/Off para encender. Seleccione un escenario en la unidad.

Paso 2: Siga las instrucciones de voz para colocar los parches en el maniquí y conectar el conector de cable.

Paso 3: Presione el botón de choque cuando se le indique, luego comience la RCP según las instrucciones.

Product Name: AED Trainer
Product No.: WL120 Practi-TRAINER Essentials



Manufactured for: WNL Products 275 Hopping Brook Road,
Holliston, MA 01746
800-884-9629
www.wnlproducts.com

Date: 2025-11-20
Doc. No.: ENNE-120D-073-011
Version: C1